



NxFl ex T9001, 1/8" \* 1/4" \* 24" L WITH MALE LUER/CAP

NxFl ex T9001, 1/4" \* 7/16" \* 6" L WITH MPC BODY-INSERT/CAP

EVA TUBING, 1/4" \* 3/8" \* 3" L (3)

FILM: NxFl ex F1000

0	17 Nov., '09	ISSUED	J. Han	CHK'D BY: Z. Wu
REV	DATE	DESCRIPTION	BY	APP'D BY:

RIM Bio, Inc	
METRIC	ENGLISH
	.X ± 0.5"
	.XX ± 0.25"
THIRD ANGLE PROJECTION	

10L Bioprocess Container

RCBB-2002Q

SHEET NO:	DATE:
1 OF 1	17 Nov., '09

PART NO:	REV:
DRAWING NO:	0
RCBB-2002Q	