



26.7"

20.0"
5.0" 6.5" 7.5"

FILM: NXFlex P1000

NXFlex T9001, 1/8" * 1/4" * 36" L WITH MALE LUER/CAP

EVA TUBING, 3/8" * 1/2" * 3" L

INJECTION PORT

EVA TUBING, 1/4" * 3/8" * 3" L

NXFlex T9001, 3/8" * 5/8" * 6" L WITH MPC BODY-INSERT/CAP

20L Bioprocess Container

RIM Bio, Inc

METRIC ENGLISH

THIRD ANGLE PROJECTION		

RCBB-2001Q

SHEET NO: 1 OF 1
DATE: 24 Nov., '09

PART NO:

DRAWING NO: RCBB-2001Q
REV: 0

ISSUED	J. Han	CHK'D BY: Z. Wu
DATE	24 Nov., '09	SCALE: N/A
DESCRIPTION		
BY	J. Han	
DATE	24 Nov., '09	
REV	0	