



FILM: NXFLEX F1000C1

TUBING, C-FLEX EQUIVALENT T9001 1/8" x 1/4" x 24" L WITH FEMALE LUER AND CAP

EVA TUBING, 1/4" x 3/8" x 3" L

10.0"

14.3"

Ø .40"

8	08 Mar., '12	CHANGES PER CUSTOMER	J. Han	SCALE: N/A	1
7	07 Mar., '12	ISSUED	J. Han	SCALE: N/A	2
REV	DATE	DESCRIPTION	BY	APP'D BY:	3

RIM Bio, Inc.		2L Bioprocess Container	
METRIC	ENGLISH		
	.X ± 0.2		
	.XX ± 0.25		
THIRD ANGLE PROJECTION			
DWG BY: J. Han	CHK'D BY: Z. Wu		
APP'D BY:			

RCBB-3003Q		RCBB-3003	
SHEET NO: 1 OF 1		DATE: 22 Mar., '10	
PART NO: RCBB-3003		DRAWING NO: RCBB-3003Q	
REV: 8			