



FILM : NxFlex F1000

EVA TUBING 1/4" x 3/8" x 3" L

INJECTION PORT

EVA TUBING, 3/8" x 1/2" x 3" L

NxFlex 19001, 1/8" x 1/4" x 36" L WITH MALE LUER/CAP

0	08 Apr., '11	ISSUED	J. Han	
REV	DATE	DESCRIPTION	BY	APP'D BY:

RIM Bio, Inc		50L Bioprocess Container	
METRIC	ENGLISH		
THIRD ANGLE PROJECTION			
DWG BY: J. Han	CHK'D BY: Z. Wu	RCBB-2006Q	
APP'D BY:	SCALE: N/A	SHEET NO: 1 OF 1	
		DATE: 08 Apr., '11	
		PART NO: RCBB-2006	
		DRAWING NO: RCBB-2006Q	
		REV: 0	