



NxFl ex T9001, 1/8" * 1/4" * 24" L WITH MALE LUER/CAP

NxFl ex T9001, 1/4" * 7/16" * 6" L WITH MPC BODY INSERT/CAP

INJECTION PORT

EVA TUBING, 1/4" * 3/8" * 3" L (3)

FILM: NxFl ex F1000

RIM Bio, Inc	
METRIC	ENGLISH
	.XX ± 0.25"
THIRD ANGLE PROJECTION	

5L Bioprocess Container

RCBB-2003Q

0	17 Nov., '09	ISSUED	J. Han	CHK'D BY: Z. Wu
REV	DATE	DESCRIPTION	BY	SCALE: N/A

SHEET NO:	DATE:
1 OF 1	17 Nov., '09

DRAWING NO:	REV:
RCBB-2003Q	0

A B C D E F 8 7 6 5 4 3 2 1